

Contents

Summary of Statistics	p. 2
Raw Data Tables	p. 2
Percent Recovery Data Tables	p. 3
Z-Score Data Tables	p. 3
Z-Score Graph	p. 4
Gaussian Normal Distribution Plots	p. 7

TERMS AND CONDITIONS OF USE

Emerald Test™ does not grant exclusive use rights with respect to this material. Emerald Test™ should be acknowledged as the source of its material. It is unlawful to falsely claim copyright or other rights in Emerald Test™ material. All rights reserved. This document or parts thereof may not be reproduced in any form, stored in any retrieval system, or transmitted in any form by any means—electronic, mechanical, photocopy, recording, or otherwise—without prior written permission of the publisher, except as provided by United States of America copyright law.



3 analytes of interest

Acceptance Criteria: 40% Recovery

2 Participating laboratories

Summary of Statistical Data

Component	Consensus Mean	Standard Deviation	Deviation %	Assigned Value	Mean Recovery %
(-)-Δ8-THC	0.6	0.06	12.65	0.51	119.37
Total CBD	82.32	7.32	9.99	73.3	112.31
Total THC	2.44	0.26	12.12	2.12	115.09

Component	# Accepted	% Accepted	# Not Accepted	% Not Accepted	# Outliers	# Participants
(-)-Δ8-THC	2	100.0	0	0.0	0	2
Total CBD	2	100.0	0	0.0	0	2
Total THC	2	100.0	0	0.0	0	2

Raw Data Tables

Analyst No.	(-)-Δ8-THC %	Total CBD %	Total THC %	Lab ID number %
3472	0.559	77.143	2.258	3472
3489	0.6495	87.4995	2.6211	3489

Percent Recovery Data

Analyst No.	(-)- Δ 8-THC	Total CBD	Total THC
3472	110.47	105.24	106.51
3489	128.36	119.37	123.64


Z-Score Data Tables


Analyst No.	(-)- Δ 8-THC	Total CBD	Total THC
3472	-1.0	-1.0	-1.0
3489	1.0	1.0	1.0




Z-Score Graphs

Total CBD Z-Score











Gaussian Normal Distribution Plots

Total CBD Distribution



Total THC Distribution



(-)- Δ 8-THC Distribution

