

Contents

Summary of Statistics	p. 2
Raw Data Tables	p. 3
Percent Recovery Data Tables	p. 4
Z-Score Data Tables	p. 5
Z-Score Graph	p. 6
Gaussian Normal Distribution Plots	p. 7

TERMS AND CONDITIONS OF USE

Emerald Test™ does not grant exclusive use rights with respect to this material. Emerald Test™ should be acknowledged as the source of its material. It is unlawful to falsely claim copyright or other rights in Emerald Test™ material. All rights reserved. This document or parts thereof may not be reproduced in any form, stored in any retrieval system, or transmitted in any form by any means—electronic, mechanical, photocopy, recording, or otherwise—without prior written permission of the publisher, except as provided by United States of America copyright law.

Summary of Statistical Data

Component	Consensus Mean	Standard Deviation	Deviation %	Assigned Value	Mean Recovery %
Water activity	0.96	0.02	2.16	0.97	99.48

Component	# Accepted	% Accepted	# Not Accepted	% Not Accepted	# Outliers	# Participants
Water activity	17	100.0	0	0.0	2	17

Raw Data Tables

Analyst No.	Water activity %	Lab ID number %
3262	0.973	3262
3276	0.995	3276
3307	0.844	3307
3365	0.9134	3365
3378	0.94	3378
3405	0.994	3405
3424	0.947	3424
3446	0.957	3446
3463	0.964	3463
3465	0.8092	3465
3475	0.96	3475
3478	0.9781	3478
3491	0.9484	3491
3499A	0.975	3499A
3499B	0.98	3499B
3514	0.9669	3514
3810	0.9764	3810



Percent Recovery Data


Analyst No.	Water activity
3262	100.31
3276	102.58
3307	87.01
3365	94.16
3378	96.91
3405	102.47
3424	97.63
3446	98.66
3463	99.38
3465	83.42
3475	98.97
3478	100.84
3491	97.77
N/A	100.52
N/A	101.03
3514	99.68
3810	100.66

Z-Score Data Tables

Analyst No.	Water activity
3262	0.41
3276	1.47
3307	Outli./NR
3365	-2.47
3378	-1.19
3405	1.42
3424	-0.85
3446	-0.37
3463	-0.03
3465	Outli./NR
3475	-0.22
3478	0.66
3491	-0.78
3499A	0.51
3499B	0.75
3514	0.11
3810	0.57



Z-Score Graphs



Gaussian Normal Distribution Plots

Water activity Distribution

